Lumbar support belt

Art. no. 44-7019 Model B711





Please read the entire instruction manual before using the product and then keep it for future reference. We accept no liability for any errors in the text or images and we reserve the right to make any necessary changes to technical data. If you have technical problems or other questions, please contact our customer service team.

Safety

- The belt should only be used by people with serious problems or injuries to the lower back, for example caused by tumours, bleeding problems or mental health conditions.
- The product must be used with caution by pregnant women and people with recent severe lumbar injuries or severe osteoporosis.
- After surgery or if you are suffering from other medical conditions, you should only use the product on a doctor's advice.
- You should not wear the product during bed rest or for one hour before and after meals to avoid any discomfort caused by long-term use or a full stomach. We recommend that you wear it for no more than two hours at a time and for a total of no more than six hours per day. If you are making long-term use of the product, you should take a break of 30 minutes between each period of use. Decide how long you are going to use the product for and how often you use it on the basis of your individual needs.

Using the product

1. Make sure that the belt is in the correct position. The support straps should be on the outside and the belt should be fastened around your stomach.
2. Hold the belt against your stomach on the left-hand side and stretch the right-hand side of the belt across your body. Attach the right-hand side of the belt to the left-hand side using the hook-and-loop fastening. **Note:** Make sure that the support straps are loose. You should tighten them once the belt feels comfortable.
3. Tighten the support straps. You should feel firm pressure on your back. Fix the support straps over your stomach with the hook-and-loop fastening. **Note:** The belt should give you support and a feeling of stability. If your breathing is restricted or you feel uncomfortable, the belt is too tight. Repeat steps two and three until the belt feels as if it is firmly fastened but is still comfortable. You must decide on this for yourself.

Cleaning and maintenance

Storage

- Do not leave the hook-and-loop fastener open when the belt is not in use to prevent it from catching on clothes and other objects.
- Do not store the belt with objects that have sharp edges or hooks as they can catch in the fibres of the belt.

Washing instructions

Note: The product should be washed separately as the hook-and-loop fastening can damage other clothes.

1. Hand wash in warm water using a neutral detergent. Rinse thoroughly in clean water to remove any remaining detergent. Do not wring the belt out. Instead, squeeze the water out carefully or use a dry towel to absorb the excess moisture.
2. Dry flat in a cool, well-ventilated place. Do not expose the product to direct sunlight. Do not iron. Correct care can extend the lifetime of the product.

Responsible disposal

The product should be disposed of in accordance with local regulations. If you are unsure of what to do, contact your local council.

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