# Bathroom Scale

Art.no 34-5062

Model RTC3010

Please read the entire instruction manual before using and save it for future use. We apologise for any text or photo errors and any changes of technical data. If you have any questions concerning technical problems please contact our Customer Service Department.

## Safety

- This product is a precision instrument and should be handled with care and protected from blows and rough treatment.
- Place the scale in the same place every time you weigh yourself.
- For accurate readings, the scale should be placed on a hard, level surface. An uneven floor beneath the scale can cause the readings to fluctuate.
- Weigh yourself without clothing at the same time of day, preferably in the morning before breakfast. Place your feet in the same place every time you weigh yourself.
- The displayed readings are rounded up/down to the nearest 100 g.
- The scale can be cleaned using a damp cloth, if and when required. Make sure that no moisture gets inside the product. Never use strong solvents or detergents.
- Keep the scale away from water and splashes; otherwise the electronic components can be damaged.

- Do not store the scale standing on its edge as the display will be on all the time and drain the battery.
- Remove the battery from the scale if it is not to be used for an extended period.
- The scale can be affected by changes in temperature or by being moved. Inaccurate readings may occur if the scale is too cold when used. Allow the scale to stand and acclimatise for approx. 24 hours before using it for the first time.

## Battery

Remove the isolating plastic tab between the battery and the battery holder before first use.

Changing the battery: Carefully press the locking tab in to release the battery cover and then remove the battery. Insert a new battery (CR2032). Replace the battery cover.

When the battery needs changing, "Lo" appears on the display.

## Use

**Warning!** Bear in mind that the glass surface might become slippery in humid atmospheres. Make sure that the surface is dry before stepping onto the scale.

1. Place the scale on a hard and flat surface (avoid placing the scale on a mat or other soft surface).
2. To activate the scale: Tap the top of the scale gently with one of your feet.\* After a few seconds "0.0" will appear on the display. The scale is now ready for use.

3. When "0.0" appears on the display, carefully mount the scale with your weight evenly distributed between both legs and between the toes and heels of your feet. Remain stationary until the display shows a steady reading for a few seconds. The display flashes when a reading has been taken.
4. The scale automatically shuts off when you step off it.

Note that the readings can be given in either kilograms (kg), pounds (lb) or stones (st). There is switch on the bottom to change the units.

\* The scale is activated by a built-in vibration sensor.

## Care and maintenance

Clean the product by wiping it with a soft, damp cloth. Never use solvents or abrasive cleaning agents.

## Disposal

Follow local ordinances when disposing of this product. If you are unsure of how to dispose of this product, please contact your local authority.

## Specifications

<b>Weight capacity</b>	150 kg / 330 lb
<b>Graduations</b>	100 g / 0.2 lb
<b>Battery</b>	1 x CR2032 (included)

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