# Kids Electric Toothbrush

Art.no 34-4128-1     Model PSA-1827B  
34-4128-2

Please read the entire instruction manual before using the product and save it for future reference. We reserve the right for any errors in text or images and any necessary changes made to technical data. If you have any questions regarding technical problems please contact Customer Services.

## Safety

- Only use the provided accessories.
- This toothbrush is only intended to be used for the purposes described in this manual.
- This toothbrush must not be immersed in water. If this should occur, thoroughly dry the batteries and the battery compartment.
- This toothbrush is not intended for children under 6 years old. Children using this toothbrush should be closely supervised.
- Keep this toothbrush away from hot surfaces.
- Do not use the toothbrush if it has been damaged.
- Do not disassemble or modify the product by yourself. Repairs should only be carried out by a qualified tradesman.
- Never let children play with the batteries.
- Never try to charge non-rechargeable batteries.
- If the batteries start to leak, take them out of the toothbrush immediately. Avoid contact with eyes and skin.

## Product description

Kids Electric Toothbrush with 2 brush heads and a 2-minute hourglass.

## Inserting the batteries

1. Remove the battery cover on the bottom by pulling it away from the handle.
2. Remove the old batteries.

3. Insert new AA/LR6 batteries. Note the polarity markings on the battery cover.
4. Replace the cover.

### Note:

Always change both batteries at the same time. Never mix old and new batteries. If the toothbrush is not to be used for an extended period, remove the batteries to avoid any possible risk of leakage.

## Use

- This electric toothbrush is intended for use by children over the age of six.
- Do not press too hard or scrub, simply let the brush do all the work.

## Changing the brush head

1. Twist the brush head anticlockwise and remove it from the handle.
2. Replace the brush head by twisting it clockwise.

## Care and maintenance

Rinse the brush head in running water.

## Disposal

Follow local ordinances when disposing of this product. If you are unsure of how to dispose of this product, please contact your local authority.

## Specifications

Battery: 2 x AA/LR6 (not included)

## GREAT BRITAIN

For consumer contact, please visit [www.clasohlson.co.uk](http://www.clasohlson.co.uk) and click on customer service.

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