Massage Unit

Art.no: 18-4009 Model: JL-3001-UK
34-6712 JL-3001

- Deep-acting massage with or without infrared heat for quick relief of muscle tension and pain.
- Variable speed control and ergonomic handle.
- 2 accessories for different types of massage supplied.

Please read the entire instruction manual before using and save it for future reference. We apologise for any text or photo errors and any changes of technical data. If you have any questions concerning technical problems please contact our Customer Service Department

Safety

Warning: Do not open the casing! Dangerous voltages are present in certain unprotected components inside the cover of the unit. Contact with such components can cause electric shock.

- The appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
 - Children should be supervised to ensure that they do not play with the appliance.
1. Always disconnect the unit after use, do not leave the unit unattended when the plug is connected to an electrical outlet.
 2. Do not use the massage unit in bathrooms or other damp places.
 3. Position the massage unit in such a way that it cannot fall into water or other liquids by accident.
 4. An electrical product connected to an electrical outlet is live even when not in use.
 5. Do not touch the massage unit with wet hands when it is connected to an electrical outlet.
 6. If the massage unit is connected to an electrical outlet and falls into water: **DISCONNECT THE PLUG FROM THE ELECTRICAL OUTLET** before touching the massage unit.
 7. Do not let children use the massage unit without supervision.
 8. Disconnect the unit from the electricity supply before cleaning.
 9. Do not use the product if the plug, cable or any other part of the unit is damaged.
 10. Use only the supplied accessories. Do not disassemble or modify the product.
 11. The massage unit is only intended for indoor use. Do not use it where flammable gases/vapours, such as aerosol propellants, oxygen gas etc, are present.
 12. Ensure that no foreign objects enter the openings on the unit to prevent electrical shocks.
 13. Disconnect the unit and let it cool down completely before storing.
 14. The product should only be repaired by qualified service personnel.

Operation

Warning:
Consult your doctor before using the unit if you are pregnant.

Do not massage the same area for more than 5 minutes and do not exceed 30 minutes in one session.

1. Ensure that the power switch is in the "0" position. Connect the unit to an electrical outlet
2. Push the power switch to desired mode:
 - Mode "M" for massage, control the speed using the knob from slow (min) to intensive (max).
 - Mode "M/H" for massage with infrared heat, control the speed using the knob from slow (min) to intensive (max).
3. Push the power switch to the "0" position when finished and disconnect the unit from the electrical outlet.

Tip: Drink one or two glasses of water before and after the massage, this helps to expel waste products from the muscles.

Care

Wipe the massage unit periodically with a lightly dampened soft cloth. Use a mild detergent; avoid solvents and aggressive or abrasive detergents.

Disposal

Follow local ordinances when disposing of this product. If you are unsure about how to dispose of this product contact your municipality.

Specifications

Model: JL-3001, JL-3001-UK
Rated voltage: 220/240 V AC, 50 Hz
Output: 35 W

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