Cordless Mini Massager

Art.no 44-4736 Model CM2210

Please read the entire instruction manual before using the product and then save it for future reference. We reserve the right for any errors in text or images and any necessary changes made to technical data. In the event of technical problems or other queries, please contact our Customer Services.

Safety

Do NOT use the product:

- On animals or children.
- On burn injuries, open cuts or scars.
- On pregnant women, recently operated or sensitive body parts such as the head and joints.
- If you have a pacemaker, implant or other physical aid, or if you are unsure whether the product is appropriate for you, you should consult a doctor before using the product.
- The product is only intended to be used for massaging a human body.
- The product must not be used by anyone (including children) suffering from physical or mental impairment.
The product should not be used by anyone that has not read the instruction manual unless they have been instructed in its use by someone who will take responsibility for their safety.
- Never let children play with the product.
- The product is only intended for personal domestic use and only in the manner described in this instruction manual.
- Never use accessories not approved by the manufacturer.
- Never attempt to modify the product in any way.
- Never use the product if it is damaged or functions abnormally.
- Keep the product away from hot surfaces.
- The product must not be immersed in water or other liquid.

Health considerations

- Contact a doctor before using this product if you are pregnant, use a pacemaker or have any of the following medical diagnoses: Diabetes, heart disease, malignant melanoma or a condition which requires substantial rest.
- If soreness, rashes or pain result after use, discontinue use immediately and contact your doctor.
- Never use the product on swollen or inflamed areas. Do not use the product on skin lesions or open wounds.
- Do not use the massager for longer than 15 minutes at a time.

Charging

Connect the included charging cable to the massager and to a USB-charger or computer.  will shine a steady red while charging is in progress and change to blue when the battery is fully charged.

Instructions for use

- Short-press  to switch the massager on. Short-press repeatedly to change the speed of the massager (4 speeds).
- Long-press  to switch the massager off.

Care and maintenance

Clean the product by wiping it with a damp cloth. Use only mild cleaning agents, never solvents or corrosive chemicals.

Responsible disposal

This symbol indicates that the product should not be disposed of with general household waste. This applies throughout the entire EU. In order to prevent any harm to the environment or health hazards caused by incorrect waste disposal, the product must be handed in for recycling so that the material can be disposed of in a responsible manner. When recycling the product, take it to your local collection facility or contact the place of purchase. They will ensure that the product is disposed of in an environmentally sound manner.

Specifications

Battery	3.7 V, 900 mAh (built-in, non-replaceable)
Power	1.5 W
Charge input	USB, 5 V DC
Operating temperature	5–40 °C
Dimensions	10 × 8 × 7.5 cm
Weight	180 g

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