Light therapy body brush

Art. no. 44-6904 Model HML-BS-WI6

Please read the entire instruction manual before use and save it for future reference. We reserve the right for any errors in text or images and for making any necessary technical changes to this document. If you should have any questions concerning technical problems please contact our Customer Services.

Safety

- The product is not intended for use as a medical treatment.
- The product must not be used by children under the age of 12.
- This appliance can be used by children aged from 13 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.
- Adequate precautions must be taken to safeguard the user against the dangers of excessive exposure to visible and infra-red radiation.
- If unexpected side effects, such as itching, occur within 48 hours of the first session of using the appliance, medical advice should be sought prior to further exposure.
- This appliance contains batteries that are only replaceable by skilled persons.
- This appliance contains batteries that are non-replaceable.
- Never use the product if it is malfunctioning or is damaged in any way. Do not use the product if you have dropped it on the floor. The appliance must not be used if the light panel or silicone surface is cracked, broken, or removed.
- Never stand or sit on the product.
- Never put the product close to a heat source, such as a cooker, woodburning stove, oven or open fire.
- Never cover the product while it is in use.
- Never try to open the housing or repair or modify the product in any way.
- **Risk for eye damage.** Do not look directly at the LEDs.
- The product is not suitable for use by anyone who has any of the following conditions or who matches any of the following descriptions:
 - Silicone allergy
 - Heart disease
 - Mental illness
 - Abnormal blood pressure
 - Sunburnt skin
 - Keloid scars
 - Oversensitivity to light
 - Discolouration and contact dermatitis
 - Infectious diseases
 - Immunological diseases
 - Diabetes
 - Melanomas on the face
 - Abnormal wound healing
 - Malignant tumours
 - People who are taking medication, cancer patients, especially with skin cancer or precancerous moles
 - People who have had botox injections in the face within the last three months
 - Pregnant; breastfeeding; and other sensitive people
 - Patients with capillaries caused by long-term use of steroid hormone medications or by liver dysfunction
- If you have permanent or temporary fillers in the area that will be covered, consult a doctor before using the product.
- If you are undergoing hospital treatment, consult a doctor before using the product.
- If you feel a burning sensation or other discomfort during the treatment, stop the treatment immediately.

Cleaning and maintenance

Rinse the brush under running water. Never use solvents, abrasive cleaners or corrosive chemicals. Be careful not to damage the light panel or silicone surface.


Simplified EU Declaration of Conformity

Clas Ohlson AB hereby declares that this type of radio equipment **HML-BS-WI6** complies with Directive 2014/53/EU. The full text of the EU declaration of conformity is available on the following website: www.clasohlson.com

Ver. 20260615

CAPERE

Product description

1. Lights
2. [O] On/Off/Light mode selector
3. LED indicator
4. USB C port
5. Charging stand

Charging

Charge the product fully before using it for the first time.

1. Connect a USB cable to the USB C port (4) on the charging stand and to a suitable USB charger.
2. Place the brush on the charging stand.
3. The LED indicator (3) shines with a red light while charging is in progress. When the battery is fully charged, the LED-indicator shines with a green light.
4. The LED indicator (3) will shine with a yellow light when the battery needs charging.

Photosensitivity test

Note You should carry out a photosensitivity test before using the product for the first time:

Switch the product on and put it on the inside of your forearm. Leave it there for 10 minutes. If any type of redness or rash appears, we advise against using the product, as your skin is probably extra sensitive. Consult a skin specialist before continuing to use the product.

Using the product

1. Hold [O] (2) for 2 seconds to start the brush. It starts with red light.
2. Press [O] repeatedly to toggle between the following modes:
 Red light (660 nm) → Blue light (450 nm) →
 Purple light (660 + 450 nm) → Purple light with vibration ↩
3. Hold [O] (2) for 2 seconds to turn the brush off.

Recommended Schedule of Exposure

Application Method: Apply the silicone brush surface directly onto the target skin area.

Max. Session Duration: Do not exceed 30 minutes per session.

Daily Frequency: This product is designed for daily care and can be used up to 2 times per day (e.g., morning and evening routines).

Responsible disposal

This symbol indicates that the product should not be disposed of with general household waste. This applies throughout the entire EU. To prevent any harm to the environment or health hazards caused by incorrect waste disposal, the product must be handed in for recycling so that the materials can be disposed of in a responsible way. To recycle the product, take it to your local recycling facility or contact the retailer. They will ensure that the product is disposed of in a way that does not harm the environment.

Specifications

Input	5 V DC/2 A
Frequency bands	110.8-141,8 kHz
Max. output power	-13.9 dBua/m @10m
Charging port	USB-C
Charging cable	USB-C – USB-C (included)
Battery	2000 mAh Li-po
Charging time	approx 4 h
Battery life	up to 2 h
IP rating	IP65
Wavelength	Red: 660 nm Blue: 450 nm Purple: 660 nm + 450 nm
Dimensions	Brush: Ø 120 mm h: 28 mm Charging stand: 86 × 85 × 73 mm
Weight	Brush: 188 g Charging stand: 57 g

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